# WHS0850 QA Freeze Report


Aptean Ltd Copyright © 2011-2025.

## **Contents**

1 WHS0850 QA Freeze Report......1

### 1 WHS0850 QA Freeze Report

This is the QA Freeze Report (WHS0850 v1.4) screen.



#### **Generating a QA Freeze Report:**

The screen displays the following input fields:

- Owner. Enter the Owner code. When selected, the owner name will be displayed next to it. An LOV is provided.
- Owner Product Code From/To: Input fields for entering a range of owner product codes. An LOV is provided.
- Reference Number From/To: Input fields for entering a range of reference numbers. An LOV is provided.
- Batch From/To: Input fields for entering a range of batch numbers. An LOV is provided.
- Manu Date From/To: Input fields for entering a range of manufacturing dates.
- Detail/Summary: A radio button group to select the report's level of detail. Options are Detail and Summary.
- Order by: A group of radio buttons to determine the report's ordering. Options are Product Code and Freeze Reason Code.

There is also a button:

#### Report

When details are entered, click the **Report** button, which will start QA Freeze Report generation. If this is a report, this will ask you how you would like to print the report through the standard screen. More details on this process can be found in the guide Printing From the WMS.

